CareFusion

# MicroDL

Operating Manual - English

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## Contents

I.	Introduction – English.....	3
	Package contents .....	4
	Configuring the MicroDL .....	5
	Configuring the MicroDL with SPCS.....	27
	Battery Management .....	28
	Calibration Check .....	29
	Cleaning the Transducer.....	30
	Consumables/ Supporting products .....	31
	Training .....	31
	Specification of the MicroDL.....	33
II.	Introduction – CareFusion USA .....	35
	Package contents .....	36
	Configuring the MicroDL.....	38
	Configuring the MicroDL with SPCS.....	61
	Battery Management .....	63
	Calibration Check .....	63
	Cleaning the Transducer.....	65
	Specification of the MicroDL.....	66
III.	Einleitung – Deutsch.....	68
	Packungsinhalt.....	69
	Konfigurierung des MicroDL .....	71
	Konfigurierung des MicroDL mit SPCS .....	96
	Batteriemanagement .....	97
	Überprüfung der Kalibrierung.....	98
	Reinigung des Messwertwandlers .....	100
	Spezifikationen des MicroDL .....	101

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## Operating manual

### **I. Introduction - English**

The MicroDL spirometer is a compact, battery operated, and fully portable data recording spirometer designed to collect spirometry data, expiratory Flow/Volume and Volume/Time curves.

The MicroDL can be configured for use with a single patient for detailed monitoring over a period or for collecting spirometry data from a number of patients.

Alternatively, the data may be uploaded to a PC running the fully featured Spirometry PC Software (SPCS). Once uploaded, the data may be viewed individually or as a trend graph and the Flow/Volume loop of any individual point examined. SPCS may also be used to perform live spirometry tests and show the expiratory Flow/Volume loop, Volume/Time curve and child incentive displays.

In single patient mode the user can perform a spirometry test. After each spirometry manoeuvre the results are displayed together with a 'traffic light' indication of their performance as an aid to asthma management. The levels that activate the red, amber and green lights, as a percentage of the predicted value, can be selected when the unit is configured. Alternatively, the activation may be based on the patient's personal best values when the unit is configured using the SPCS software.

In multiple patient mode spirometry results can be collected from up to 200 patients. In this mode the 'traffic light' activation is based upon the patient's predicted values.

The spirometer uses the CareFusion Digital Volume Transducer, an extremely stable form of volume transducer, which measures expired air directly at B.T.P.S (Body Temperature and Pressure with Saturated water vapour) thus avoiding the inaccuracies of temperature corrections. This transducer is insensitive to the effects of condensation and temperature and avoids the need for individual calibration prior to performing a test.

## **Package contents**

The MicroDL is packaged in a convenient carrying case and comes complete with the following items (Fig.1):

1. MicroDL microcomputer unit
2. CareFusion Digital Volume Transducer
3. PP3 alkaline battery
4. Serial Cable together with disposable cardboard mouthpieces.


Fig. 1

## Configuring the MicroDL


The MicroDL must be configured for single patient (asthma management) or multiple patient use. Multiple patient use is configured by hand on the instrument itself and single patient use may be configured either by hand or more conveniently by the SPCS software provided (see SPCS installation and operation).

When configuring for single patient mode a patient identifier must be entered. This number will be used to identify data when printing or uploading results. The age, height and sex of the patient together with any ethnic correction must also be entered. This data is used to calculate the predicted values for FEV1, FVC, PEF, and FEV1/FVC. The ethnic correction reduces the predicted values of the volume measurements by the percentage entered. Finally, the levels at which the amber and red warning lights are activated are entered.

These levels are expressed as a percentage of the predicted values when the unit is configured by hand but the patient's personal best values may be used when the unit is configured using SPCS.

When using the MicroDL in multiple patient mode the patient details are entered each time the unit is turned on. Therefore, when configuring the unit for multiple patient use, only the warning light levels are entered.


To configure the unit, depress the off key, turn the unit on and release the off key. If there is data currently stored in the unit then the number of records displayed:

Configuring the unit will delete these records so if the stored data is required turn the unit off and print or upload the data using SPCS (see Printing Stored Records and Uploading Results).


To continue with configuration and delete the data, press 'enter' and the following will be displayed:


The 1 displayed indicates that the unit is set to single patient mode.

By using the up or down key, as indicated by the flashing arrows, multiple patient mode may be selected when the display will change to:


Once the required configuration is obtained, press 'enter'.

If multiple patient mode has been selected the amber light will flash and the display will change to:


Use the up and down keys to select the level of results below which will produce an amber warning light when the user performs a spirometry manoeuvre. This level is expressed as a percentage of the predicted value. Please note that the default values displayed may vary from those shown in this manual as the MicroDL will automatically update the default values to the last actual values used.

After pressing 'enter', the red light will flash and the display will change to:

Use the up and down keys to select the required percentage of predicted value which will produce a red warning light when the user performs a spirometry manoeuvre less than the set value.

When the required level has been selected press 'enter' and the display will change to:

The unit is now configured for multiple patient use and will automatically turn off.

The most convenient way to configure the unit for a single patient is to use the SPCS software provided. SPCS allows the user to configure several options that are not available when the unit is configured by hand.

These include:


Personal best values can be entered as a basis for the warning light operation.

A simple operating procedure, for use by the patient, can be printed.


The results to be displayed can be selected.

However, the unit may be configured by hand in the following manner.

Turn the unit on with the 'off' key depressed as described above and select single patient mode by pressing 'enter' when the display shows:


The display will then change to:


Use the up and down keys to select a patient number.


This number can be set using the 'up' and 'down' keys anywhere in the range 1 to 1999 and will be used to identify data when printing or uploading results. When the required number is displayed press 'enter' and the display will change to:

The flashing symbol represents a male patient and using the up or down key will change the display to indicate a female patient:

When the correct selection has been made press 'enter' and the display will change to:


Adjust the patients age using the up and down keys and press 'enter'. The display will change to:

Enter the patients height, in centimetres, using the up and down keys and press 'enter'.

The display will change to:


Adjust the ethnic correction using the up and down keys. The display shows the percentage reduction of volumetric predicted values to be applied for the patient. Press 'enter' when the required correction has been selected. The levels for the amber and red warning lights should now be set as described above in the description of multiple patient configuration.


## Spirometry Operation (Single patient use)

When performing a spirometry test the user should be instructed to insert a disposable or clean reusable, mouthpiece into the mouthpiece holder of the spirometer.

The unit is then turned on, the current time will be displayed briefly, and the three lights will be illuminated before the display changes to:

The patient should then inhale as deeply as possible, seal their lips around the mouthpiece and exhale as hard and as fast as possible until no more air can be exhaled. The spirometry results are then displayed:

The first result shown is Peak Expiratory Flow (PEF). The arrows indicate that the up and down keys may be used to scroll through the other results of the test. These include FEV1, FVC, FEV1/FVC (shown as FEV1 FVC), followed by the percentage of predicted value for all of the above.

Note: when the unit is configured with SPCS software the displayed results can be limited to only those that the patient can usefully interpret, for example, the unit can be configured to display only PEF and the percent predicted PEF. However, regardless of the

measurements displayed, all measurements (including expiratory Flow/Volume and Volume/Time data) will be recorded.


In addition, the MicroDL can be configured with PEF displayed in units of Litres/min or Litres/sec.

As well as the spirometry results, warning lights will be displayed. The lights are based on a percentage of the predicted value or personal best values.

Note: personal best values can only be used when the unit is configured using SPCS software.

The percentage level below which the warning lights illuminate can be configured using SPCS software or when configuring the unit by hand. The warnings lights will vary according to the spirometry parameter displayed. However, if a lower level was present on one of the other parameters then the lower level light will flash until the display is scrolled to that parameter.

For example, if the unit is configured to display PEF, FEV<sub>1</sub>, and FVC, then the patient will see the following displays when scrolling through the results after a manoeuvre with good peak flow but a poor FEV<sub>1</sub> and very poor FVC:

After completing a spirometry manoeuvre the user may press 'enter' to repeat or off to finish testing. After completing a sequence of tests, the best result will be stored. The definition of best is selected using SPCS and may be best FEV1, best FVC, best PEF, or best FEV1 +FVC. If the unit is turned off after performing a manoeuvre, and then used again within five minutes, then only the best result within the period will be saved.

## **Spirometry Operation (Multiple patient use)**


When performing a spirometry test the user should be instructed to insert a disposable or clean reusable, mouthpiece into the mouthpiece holder of the spirometer.

The unit is then turned on, the current time will be displayed briefly, and the three lights will be illuminated before the display changes to:


Use the up and down keys to select a patient number. This number can be set, using the up and down keys, anywhere in the range 1 to 1999 and will be used to identify data when printing or uploading results.


When the required number is displayed press 'enter' and the display will change to:

The flashing symbol represents a male patient and using the up key will change the display to indicate a female patient:


Pressing the up key again will remove the patient symbol:

This indicates that predicted values are not required and the unit will not request the patient's details. Please note that if this option is selected then the warning lights will not be operational.

If male or female predicted values have been selected then the display will change to:


Select the patients age using the up and down keys and press 'enter'. The display will change to:


Enter the patients height, in centimetres, using the up and down keys and press 'enter'. The display will change to:

Adjust the ethnic correction using the up and down keys. The display shows the percentage drop of volumetric predicted values to be applied for the patient. Press 'enter' when the required correction has been selected.

Note: the initial patient number displayed and the default patient details will be the same as the last one used. This is to facilitate quick entry of details when performing a repeat test on the same patient. If a new patient is to be tested then select a different patient number.


After the patient details have been entered, the display will change to:

Insert a disposable or clean reusable, mouthpiece into the mouthpiece holder of the spirometer.

Instruct the patient to inhale as deeply as possible, seal their lips around the mouthpiece and exhale as hard and as fast as possible until no more air can be exhaled.

The spirometry results are then displayed:

The first result shown is Peak Expiratory Flow (PEF). The arrows indicate that the up and down keys may be used to scroll through the other results of the test.


These include FEV1, FVC, FEV1/FVC (shown as FEV1 FVC), followed by the percentage of predicted value for all of the above.

Note: if the unit has been previously configured with SPCS, the displayed results will be limited to only those previously set. However, regardless of the measurements displayed, all measurements (including expiratory Flow/Volume and Volume/Time data) will be recorded.

As well as the spirometry results, warning lights will also be displayed as described in the section 'Spirometry Operation (Single patient use)'. However, in the case of multiple patient use, the light activation is based upon the percentage of the predicted value as personal best values cannot be entered.

After completing a spirometry manoeuvre the user may press 'enter' to repeat or 'off' to finish testing. After completing a sequence of tests, the best result will be stored. The definition of best is selected using SPCS and may be best FEV1, best FVC, best PEF, or best FEV1+FVC. When off is pressed, the user has the option to obtain a printout of the results.

The display will change to:

If a printout is not required, press 'off' again and the unit stores the results, date and time stamped, and turns off. If a printout is required then connect the serial to parallel converter to the parallel input socket on the printer. Connect the jack plug to the socket on the side of the MicroDL and turn the printer on.

NOTE: Keep the printer out of reach of the patient when connected to the MicroDL.

**Press 'enter' on the MicroDL to obtain the following report:**

MicroDL FULL REPORT


Patient Name: \_\_\_\_\_

Test: 1      Date: 12/11/97      Time: 0802

ID: 6      Sex: Male      Age: 19

Height: 190 cm      Ethnic Correction: 0%

Parameter	Result	%Pred	Pred	Min	Max	Units
FEV1	3.71	74	4.95	4.11	5.79	L
FVC	4.61	77	5.95	4.95	6.95	L
PEF	9.67	90	10.74	8.75	12.73	L/S
FEV1%FVC	80	97	82	70	94	%
F50	4.15	68	6.07	3.90	8.24	L/S
F25	1.36	45	2.96	1.68	4.24	L/S
MEF	5.45	102	5.31	3.60	7.20	L/S
FET	4.66					S

Normal Values: ECCS (adult); Zapletal, Solymar, Cogswell (child)  
Results at BTPS

**Technician:** \_\_\_\_\_ **Physician:** \_\_\_\_\_

The first line gives the MicroDL version number and type of report (a summary report of all data stored is also available – see Printing Stored Results).

Subsequent lines give the test number (the tests are numbered sequentially as they are stored), date and time of test, and patient details.

This is followed by the test results presented with the percentage of the predicted value, the predicted value, the predicted minimum and maximum, and the unit of measurement.


The expiratory Flow/Volume and Volume/Time curves are also shown overlaid with the predicted curves and FEV1 marker on the Flow/Volume curve.

To obtain subsequent copies of the report press 'enter'. When all the required printouts have been obtained press 'off' and the unit will turn off.

If the unit is turned off after performing a manoeuvre, and then used again within five minutes, then only the best result within the period will be saved.

## Printing Stored Records

There are five printing options for the stored results. To assist the selection of the required option a system of icons will appear on the display when the unit is set to printing mode:

If one page appears, then a single record has been selected.

If two pages appear, then only records since the last printout will be printed.

If three pages appear, then all stored records will be printed.

If the front page shows only lines of text then a summary report has been selected.

If the Flow/Volume symbol appears, then a spirometry report has been selected.


Depress the up key, turn the unit on and release the up key after a few seconds. The display will indicate the type of report to be printed and the number of pages that will be used.

Use the up key to select the type of printout required.

The following examples assume that the unit is configured in single patient mode.


### **Option 1 – summary report of all records**

After turning the MicroDL on with the up key depressed the following will be displayed:

This display is indicating that a summary report of all stored results has been selected and that the printout will use two sheets of paper.

If this option is required then press 'enter' and the display will change to:

**The following printout of spirometry results will be obtained:**

MicroDL      SUMMARY REPORT						
Test	ID	Date	Time	FEV1	FVC	PEF
-----						
1	6	10/09/97	10:56	2.97(86%)	3.25(78%)	7.49(87%)
2	6	11/09/97	11:32	2.74(79%)	3.12(75%)	6.45(75%)
3	6	12/09/97	09:58	3.14(91%)	3.53(85%)	8.37(97%)
✓	✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓	✓

Values in brackets are percentage of best value.

**Technician:** \_\_\_\_\_ **Physician:** \_\_\_\_\_

The summary report gives the date, time, FEV1, FVC, and PEF together with the percentage of best value in brackets. The tests are sequentially numbered and this can be used to select a particular test for more detailed investigation by selecting a full report using printing option 2. At any time during printing pressing 'off' will terminate the printing. Should the printer fail to operate correctly then turn off both the printer and the MicroDL unit, check all connections, and try again.


**Option 2 – full spirometry report of one record**

After turning the MicroDL on with the up key depressed, use the up key again to obtain printing option 2. The display will show:


This option is used to select a full report, including Flow/Volume and Volume/Time curves, on one stored result. The '1' indicates that one sheet of paper will be used.

Press 'enter' if this option is required and the display will change to:

The tests are numbered sequentially as they are stored and the number displayed represents the last test recorded.

Use the down key to select the required test and then press 'enter'. The display will change to:

**The following printout of spirometry results will be obtained:**

MicroDL FULL REPORT

Patient Name: \_\_\_\_\_


Test: 1 Date: 12/11/97 Time: 0802

ID: 6 Sex: Male Age: 19

Height: 190 cm Ethnic Correction: 0%

Parameter	Result	%Best	Best	%Pred	Pred	Min	Max	Units
FEV1	3.71	78	4.76	74	4.95	4.11	5.79	L
FVC	4.61	82	5.62	77	5.95	4.95	6.95	L
PEF	9.67	92	10.65	90	10.74	8.75	12.73	L/S
FEV1%FVC	80	97	82	97	82	70	94	%
F50	4.15			68	6.07	3.90	8.24	L/S
F25	1.36			45	2.96	1.68	4.24	L/S
MEF	5.45			102	5.31	3.60	7.20	L/S
FET	4.66							S

Normal Values: ECCS (adult); Zapletal,Solymar,Cogswell (child)  
 Results at BTPS

**Technician:**\_\_\_\_\_ **Physician:**\_\_\_\_\_

The first line gives the MicroDL version number and type of report. Subsequent lines give the test number, date and time of test, and patient details. This is followed by the test results presented with the percentage of personal best value, the personal best value, the percentage of predicted value, the predicted value, the predicted minimum and maximum, and the unit of measurement.


The expiratory Flow/Volume and Volume/Time curves are also shown overlaid with the predicted curves and FEV1 marker on the

Flow/Volume curve. Please note that personal best values cannot be entered in multiple patient mode.

To obtain subsequent copies of the report press 'enter'. When all the required printouts have been obtained press 'off' and the unit will turn off.


### **Option 3 – full spirometry report of all records**

After turning the MicroDL on with the up key depressed, use the up key twice to obtain printing option 3. The display will show:

This option is used to print full spirometry reports, including Flow/Volume and Volume/Time curves, on all the stored results.

The '65' indicates that sixty-five sheets of paper will be used, as there are currently sixty-five tests stored. Press 'enter' if this option is required. The display will show:

The printouts obtained will be identical to the single report produced with option 2.

At any time during printing pressing 'off' will terminate the printing. Should the printer fail to operate correctly then turn off both the printer and the MicroDL unit, check all connections, and try again.

#### **Option 4 – summary report of all records since last printout**


After turning the MicroDL on with the up key depressed, use the up key three times to obtain printing option 4.

The display will show the following:

This option is used to print a summary report of all the results that have been recorded since the last summary report was printed.


Press 'enter' if this option is required. The display will show:

The printouts obtained will be identical to the summary report produced with option 1, but with fewer entries if a summary report has been previously obtained.


#### **Option 5 – full spirometry report of all records since last printout.**

After turning the MicroDL on with the up key depressed, use the up key four times to obtain printing option 4. The display will show:

This option is used to print full spirometry reports of all the results that have been recorded since the last time that printing option 3 was used.

Press 'enter' if this option is required. The display will show:

The printouts obtained will be identical to the single report produced with option 3.

## PC Based Spirometry

The MicroDL can be connected to a PC and interfaced to CareFusion's spirometry data management software SPCS. SPCS enables you to record, manage, and display spirometry information. Its main features are:

A fully featured spirometry database.

An interface, which enables the unit to be controlled from the PC and provides live display of Flow/Volume, and Volume/Time graphs together with a child incentive display.


Reporting facility including a full set of parameters, predicted values, graphs and interpretations.

Trend analysis of spirometry data.

## Configuring the MicroDL with SPCS

The MicroDL can be set up using SPCS. First, ensure that the MicroDL is connected to the PC using the cable supplied. Connect the nine way 'D' connector to the serial port of the PC and the 3 pole jack plug to the socket in the right hand side of the MicroDL.

To communicate with the PC the MicroDL must be turned on in PC mode. Depress the 'enter' key, turn the unit on, and release the 'enter' key when the display shows:

Start SPCS on the PC. Ensure that SPCS is set up to use the COM port on the computer that the MicroDL is connected to (refer to online SPCS help). A patient record should then be set up containing the details to be downloaded to the MicroDL (or an existing patient selected). These details can then be used to initialise the unit. Please refer to the online help included with SPCS for the exact details for performing these operations.

NOTE: Keep the PC out of reach of the patient when connected to the MicroDL. The MicroDL should only be connected to a computer that is manufactured in accordance with EN60950 1992/1993 'Safety of Information Technology Equipment including Electrical Business Equipment'.

## Uploading Results

To upload the results stored in the MicroDL, connect the unit to the PC as described above and start SPCS. When the initial screen is displayed, click on the 'upload' button. The data contained in the

unit will then be uploaded into the SPCS database (refer to online help for further information).

Please note that the patient details for any data uploaded into SPCS from a unit that has been set-up by hand, should be entered into the SPCS database. This will allow the uploaded spirometry records to be displayed and trended.

The patient IDs used in the unit should match those in the SPCS database.

NOTE: Keep the PC out of reach of the patient when connected to the MicroDL.

To connect the unit to the PC; connect the cable to one of the PC's serial communications ports and turn the unit on with the 'enter' key held down. The unit display should show 'PC' flashing. Now run SPCS on your PC (for more information on SPCS, refer to online help).

## **Battery Management**

There are two sources of internal battery power in the MicroDL.

The recorded data is stored in a battery-backed memory, supplied by a Lithium power source with an expected life of at least 10 years. This battery also supplies the timekeeping clock.

The microcomputer circuitry is supplied by a 9 volt, alkaline PP3 battery. The battery life is approximately 24 hours of continuous use.

The battery voltage is monitored continuously and the following warning message is displayed when the battery starts to become exhausted:

The warning will flash for a few seconds and the unit will turn itself off. Replace the battery by sliding off the battery cover at the back of the unit, unclip the exhausted battery and replace with a new battery of the correct type.

Note: if the battery is not to be used for long periods, the battery should be removed to prevent damage to the instrument by possible leakage.

## Calibration Check

The MicroDL is calibrated to read in litres at body temperature, barometric pressure, and saturated with water vapour (BTPS) for FEV1 and FVC. The units of PEF are litres per second BTPS or litres per minute BTPS depending upon the configuration in the PC software.

The calibration should remain stable indefinitely, unless the transducer is physically damaged, and the unit should not require re-calibration.

However, as a check on the correct functioning of the unit we do recommend that the calibration is checked periodically.

To check the calibration of the unit a large syringe should be used.

This should have a 30mm outlet and be of at least three litres in capacity and preferably six litres.

A 3 Litre Syringe is available from CareFusion as Cat. No. 36-SM2125. The MicroDL should be connected to this syringe with the

minimum of adapters and the syringe volume should be injected into the transducer evenly, without pausing.

If, after checking, the unit appears to require re-calibration then please refer to your supplier for service.

## **Looking after your MicroDL Spirometer**

Please observe the following precautions:

- \* Avoid exposing the MicroDL to direct sunlight.
- \* Avoid operating the spirometer in dusty conditions or near to heating appliances or radiators.
- \* Do not keep the spirometer in a damp place or expose it to extremes of temperature.
- \* Do not direct the transducer holder towards a strong light source whilst operating the spirometer.

## **Cleaning the Transducer**

The transducer requires no routine maintenance or servicing.

However, if you wish to sterilise or clean the transducer it may be removed by means of the following procedure:

1. Remove the transducer by gently pulling from the main body with a twisting action.
2. The transducer may now be immersed in warm soapy water for routine cleaning or immersed in cold sterilising solutions e.g. Perasafe, for a period not exceeding 10 minutes (Alcohol and chloride solutions should be avoided). After cleaning/sterilising, the transducer should be rinsed in distilled water and dried. Perasafe is available from CareFusion in convenient 81g plastic bottles, Cat No. 36-SSC5000A.
3. Re-assemble the transducer into the MicroDL.

## Consumables / Supporting Products

Cat. No.	Description
36-PSA1000	Adult Disposable Mouthpieces (500 per box)
36-SST1000	One-way Safety Mouthpieces (500 per box)
36-PSA1200	Paediatric Disposable Mouthpieces (250 per box)
36-PSA1100	Paediatric Adaptor
36-MGF1025	MicroGard Pulmonary Filters (25 per box)
36-MGF1100	MicroGard Pulmonary Filters (100 per box)
36-SSC5000A	PeraSafe Sterilising Powder 81g (to make up 5 litres of solution)
36-VOL2104	Nose Clips (pack of 5)
36-BAT1002	Alkaline PP3 Battery
36-SM2125	3 Litre Calibration Syringe
36-SPC1000	Spirometry PC Software

## Caution

- Pulmonary filters are single patient use. If used on more than one patient there is a risk of cross-infection. Repeat use may increase air resistance and lead to an incorrect measurement.”
- Mouthpieces are single patient use. If used on more than one patient there is a risk of cross-infection. Repeat use may degrade materials and lead to an incorrect measurement.”

## Training

CareFusion run spirometry training courses throughout the UK, sessions include:

- \* Anatomy and Physiology of Respiratory System

- \* Lung Volume Terminology
- \* Objectives of Spirometry
- \* Setting up the test & performing a good manoeuvre
- \* Errors of technique
- \* Interpretation of results
- \* Using Spirometry as an aid to diagnosing airways disease
- \* Choosing a Spirometer
- \* Case Studies

## Symbols

Type B device

In accordance with Directive 93/42/EEC

Disposal in compliance with WEEE

## Environment

This instrument complies with directive EN60601-1-2 electromagnetic compatibility but can be affected by cellular phones and by electromagnetic interference exceeding levels specified in EN 50082-1:1992

## Servicing

A full service manual including circuit diagrams and parts list is available upon request.

## Specification of the MicroDL

### Displayed Measurements:

Forced Vital Capacity (FVC)  
Forced Expired Volume in 1 second (FEV1)  
Peak Expiratory Flow Rate (PEF),  
Forced Expiratory Ratio (FEV1 FVC)

### Recorded Measurements:

Forced Vital Capacity (FVC)  
Forced Expired Volume in 1 second (FEV1)  
Peak Expiratory Flow Rate (PEF),  
Forced Expiratory Ratio (FEV1%FVC)  
Flow at 50% of volume remaining (F50)  
Flow at 25% of volume remaining (F25)  
Mid Expiratory Flow (MEF)  
Forced Expiratory Time (FET)

**Storage Capacity:** Typically > 1 year of data recording

**Display:** Custom LCD

**Transducer Type:** CareFusion Uni-Directional Digital  
Volume

**Accuracy:** To ATS recommendations  
Standardisation of spirometry

**Power supply:** 9 volt PP3 Alkaline cell

**Backup battery life:** >10 years

**Main Battery Life:** Approximately 24 hours of continuous  
use

**Dimensions:** 170 X 70 X 60mm

**Weight:** 190g

**Operating Temperature:** 0 to +40° C

**Operating Humidity:** 30% to 90% RH

**Storage Temperature:** -20 to + 70° C

**Storage Humidity:** 10% to 90% RH

# MicroDL

## Operating manual

Federal (USA) law restricts this device to sale by or on the order of a physician or licensed practitioner.

CareFusion,  
22745 Savi Ranch Parkway,  
Yorba Linda  
CA 92887-4668  
USA

For US Customer Care:  
Toll Free 1-800-231-2466, phone 1-714-283-2228  
[www.carefusion.com](http://www.carefusion.com)

## II. Introduction – CareFusion USA

The MicroDL spirometer is a compact, battery operated, and fully portable data recording spirometer designed to collect spirometry data, expiratory Flow/Volume and Volume/Time curves, and symptom scores.

The MicroDL can be configured for use with a single patient for detailed monitoring over a period or for collecting spirometry data from a number of patients.

Alternatively, the data may be uploaded to a PC running the fully featured Spirometry PC Software (SPCS). Once uploaded, the data may be viewed individually or as a trend graph with symptom scores (in single patient mode) and the Flow/Volume loop of any individual point examined. SPCS may also be used to perform live spirometry tests and show the expiratory Flow/Volume loop, Volume/Time curve and child incentive displays.

In single patient mode the user can perform a spirometry test and enter symptom scores or event markers. After each spirometry manoeuvre the results are displayed together with a 'traffic light' indication of their performance as an aid to asthma management. The levels that activate the red, amber and green lights, as a percentage of the predicted value, can be selected when the unit is configured. Alternatively, the activation may be based on the patient's personal best values when the unit is configured using the SPCS software.


In multiple patient mode spirometry results can be collected from up to 200 patients. In this mode the 'traffic light' activation is based upon the patient's predicted values and symptom scores are not used.

The spirometer uses the CareFusion Digital Volume Transducer, an extremely stable form of volume transducer, which measures expired air directly at B.T.P.S (Body Temperature and Pressure with Saturated water vapour) thus avoiding the inaccuracies of temperature corrections. This transducer is insensitive to the effects of condensation and temperature and avoids the need for individual calibration prior to performing a test.

## **Package Contents**

The MicroDL is packaged in a convenient carrying case and comes complete with the following items (Fig.1):

1. MicroDL microcomputer unit
2. Digital Volume Transducer
3. PP3 alkaline battery
4. Serial Cable together with disposable cardboard mouthpieces.

## Configuring the MicroDL


The MicroDL must be configured for single patient (asthma management) or multiple patient use. Multiple patient use is configured by hand on the instrument itself and single patient use may be configured either by hand or more conveniently by the SPCS software provided (see SPCS installation and operation).

When configuring for single patient mode a patient identifier must be entered. This number will be used to identify data when printing or uploading results. The age, height and sex of the patient together with any ethnic correction must also be entered. This data is used to calculate the predicted values for FEV1, FVC, PEF, and FEV1/FVC. The ethnic correction reduces the predicted values of the volume measurements by the percentage entered. Finally, the levels at which the amber and red warning lights are activated are entered.

These levels are expressed as a percentage of the predicted values when the unit is configured by hand but the patient's personal best values may be used when the unit is configured using SPCS.


When using the MicroDL in multiple patient mode the patient details are entered each time the unit is turned on. Therefore, when configuring the unit for multiple patient use, only the warning light levels are entered.

To configure the unit, depress the off key, turn the unit on and release the off key. If there is data currently stored in the unit then the number of records displayed:


Configuring the unit will delete these records so if the stored data is required turn the unit off and print or upload the data using SPCS (see Printing Stored Records and Uploading Results).

To continue with configuration and delete the data, press 'enter' and the following will be displayed:


The 1 displayed indicates that the unit is set to single patient mode.

By using the up or down key, as indicated by the flashing arrows, multiple patient modes may be selected when the display will change to:

Once the required configuration is obtained, press 'enter'.


If multiple patient modes has been selected the amber light will flash and the display will change to:

Use the up and down keys to select the level of results below which will produce an amber warning light when the user performs a spirometry manoeuvre. This level is expressed as a percentage of


the predicted value. Please note that the default values displayed may vary from those shown in this manual as the MicroDL will automatically update the default values to the last actual values used.

After pressing 'enter', the red light will flash and the display will change to:

Use the up and down keys to select the required percentage of predicted value which will produce a red warning light when the user performs a spirometry manoeuvre less than the set value.

When the required level has been selected press 'enter' and the display will change to:

The unit is now configured for multiple patient use and will automatically turn off.

The most convenient way to configure the unit for a single patient is to use the SPCS software provided. SPCS allows the user to configure several options that are not available when the unit is configured by hand.

These include:

Personal best values can be entered as a basis for the warning light operation.


The unit may be configured to record symptom scores and events.

A simple operating procedure, for use by the patient, can be printed.


The results to be displayed can be selected.

However, the unit may be configured by hand in the following manner.

Turn the unit on with the 'off' key depressed as described above and select single patient mode by pressing 'enter' when the display shows:


The display will then change to:


Use the up and down keys to select a patient number.


This number can be set using the 'up' and 'down' keys anywhere in the range 1 to 1999 and will be used to identify data when printing or uploading results. When the required number is displayed press 'enter' and the display will change to:

The flashing symbol represents a male patient and using the up or down key will change the display to indicate a female patient:

When the correct selection has been made press 'enter' and the display will change to:


Adjust the patient's age using the up and down keys and press 'enter'. The display will change to:

Enter the patient's height, in centimetres, using the up and down keys and press 'enter'.

The display will change to:


Adjust the ethnic correction using the up and down keys. The display shows the percentage reduction of volumetric predicted values to be applied for the patient. Press 'enter' when the required

correction has been selected. The levels for the amber and red warning lights should now be set as described above in the description of multiple patient configuration.

## Spirometry Operation (Single patient use)

When performing a spirometry test the user should be instructed to insert a disposable or clean reusable, mouthpiece into the mouthpiece holder of the spirometer.

The unit is then turned on, the current time will be displayed briefly, and the three lights will be illuminated before the display changes to:

The patient should then inhale as deeply as possible, seal their lips around the mouthpiece and exhale as hard and as fast as possible until no more air can be exhaled. The spirometry results are then displayed:

The first result shown is Peak Expiratory Flow (PEF). The arrows indicate that the up and down keys may be used to scroll through the other results of the test. These include FEV1, FVC, FEV1/FVC (shown as FEV1 FVC), followed by the percentage of predicted value for all of the above.

Note: when the unit is configured with SPCS software the displayed results can be limited to only those that the patient can usefully interpret, for example, the unit can be configured to display only PEF and the percent predicted PEF. However, regardless of the measurements displayed, all measurements (including expiratory Flow/Volume and Volume/Time data) will be recorded.


In addition, the MicroDL can be configured with PEF displayed in units of Litres/min or Litres/sec.

As well as the spirometry results, warning lights will be displayed. The lights are based on a percentage of the predicted value or personal best values.

Note: personal best values can only be used when the unit is configured using SPCS software.

The percentage level below which the warning lights illuminate can be configured using SPCS software or when configuring the unit by hand. The warnings lights will vary according to the spirometry parameter displayed. However, if a lower level was present on one of the other parameters then the lower level light will flash until the display is scrolled to that parameter.

For example, if the unit is configured to display PEF, FEV1, and FVC, then the patient will see the following displays when scrolling through the results after a manoeuvre with good peak flow but a poor FEV1 and very poor FVC:


After completing a spirometry manoeuvre the user may press 'enter' to repeat or off to finish testing. After completing a sequence of tests, the best result will be stored. The definition of best is selected using SPCS and may be best FEV1, best FVC, best PEF, or best FEV1 +FVC. If the unit is turned off after performing a manoeuvre, and then used again within five minutes, then only the best result within the period will be saved.

## Symptom Score Entry (Single patient use only)


SPCS can be used to configure the unit to record symptom scores and events. Symptom scores have a magnitude range of 1 to 10 and event markers have no magnitude. A total of four symptom scores and event markers may be recorded. The scores may be uploaded to SPCS or printed with the spirometry results.

The symptom score or event markers are recorded by pressing 'enter' when the blow legend is displayed:


The display will change to:


Note that pressing 'enter' when the E0 legend is displayed does not record an event, but returns the display to the previous spirometry screen.

Use the up and down keys to scroll to the required event number (E1 to E4) e.g:

Press 'enter' when the required event is displayed. If the selected event had no magnitude assigned to it, the event will be recorded and the display will change to show E0 again. If another event is to be recorded, then repeat the procedure above. Press 'enter' when E0 is displayed to finish recording events.

If a selected event has magnitude assigned to it then the display will change to:

The number before the decimal point (2) refers to the event number and the number on the right hand side is the response.


Use the up and down keys to adjust the response to the required value. When 'enter' is pressed the magnitude of the event will be recorded and display will show E0 again.

After all the required events have been entered the unit may be turned off or a spirometry manoeuvre may be performed.

## **Spirometry Operation (Multiple patient use)**


When performing a spirometry test the user should be instructed to insert a disposable or clean reusable, mouthpiece into the mouthpiece holder of the spirometer.

The unit is then turned on, the current time will be displayed briefly, and the three lights will be illuminated before the display changes to:

Use the up and down keys to select a patient number. This number can be set, using the up and down keys, anywhere in the range 1 to 1999 and will be used to identify data when printing or uploading results.


When the required number is displayed press 'enter' and the display will change to:

The flashing symbol represents a male patient and using the up key will change the display to indicate a female patient:


Pressing the up key again will remove the patient symbol:

This indicates that predicted values are not required and the unit will not request the patient's details. Please note that if this option is selected then the warning lights will not be operational.

If male or female predicted values have been selected then the display will change to:


Select the patients age using the up and down keys and press 'enter'. The display will change to:

Enter the patient's height, in centimetres, using the up and down keys and press 'enter'.


The display will change to:

Adjust the ethnic correction using the up and down keys. The display shows the percentage drop of volumetric predicted values to be applied for the patient. Press 'enter' when the required correction has been selected.

Note: the initial patient number displayed and the default patient details will be the same as the last one used. This is to facilitate quick entry of details when performing a repeat test on the same patient. If a new patient is to be tested then select a different patient number.

After the patient details have been entered, the display will change to:

Insert a disposable or clean reusable, mouthpiece into the mouthpiece holder of the spirometer.

Instruct the patient to inhale as deeply as possible, seal their lips around the mouthpiece and exhale as hard and as fast as possible until no more air can be exhaled.

The spirometry results are then displayed:

The first result shown is Peak Expiratory Flow (PEF). The arrows indicate that the up and down keys may be used to scroll through the other results of the test.


These include FEV1, FVC, FEV1/FVC (shown as FEV1 FVC), followed by the percentage of predicted value for all of the above.

Note: if the unit has been previously configured with SPCS, the displayed results will be limited to only those previously set. However, regardless of the measurements displayed, all measurements (including expiratory Flow/Volume and Volume/Time data) will be recorded.

As well as the spirometry results, warning lights will also be displayed as described in the section 'Spirometry Operation (Single patient use)'. However, in the case of multiple patient use, the light activation is based upon the percentage of the predicted value as personal best values cannot be entered.

After completing a spirometry manoeuvre the user may press 'enter' to repeat or 'off' to finish testing. After completing a sequence of tests, the best result will be stored. The definition of best is selected using SPCS and may be best FEV1, best FVC, best PEF, or best FEV1+FVC. When off is pressed, the user has the option to obtain a printout of the results.

The display will change to:

If a printout is not required, press 'off' again and the unit stores the results, date and time stamped, and turns off. If a printout is required then connect the serial to parallel converter to the parallel input socket on the printer. Connect the jack plug to the socket on the side of the MicroDL and turn the printer on.

NOTE: Keep the printer out of reach of the patient when connected to the MicroDL.

**Press 'enter' on the MicroDL to obtain the following report:**

MicroDL FULL REPORT


Patient Name: \_\_\_\_\_

Test: 1      Date: 12/11/97      Time: 0802

ID: 6      Sex: Male      Age: 19

Height: 190 cm      Ethnic Correction: 0%

Parameter	Result	%Pred	Pred	Min	Max	Units
FEV1	3.71	74	4.95	4.11	5.79	L
FVC	4.61	77	5.95	4.95	6.95	L
PEF	9.67	90	10.74	8.75	12.73	L/S
FEV1%FVC	80	97	82	70	94	%
F50	4.15	68	6.07	3.90	8.24	L/S
F25	1.36	45	2.96	1.68	4.24	L/S
MEF	5.45	102	5.31	3.60	7.20	L/S
FET	4.66					S

Normal Values: ECCS (adult); Zapletal,Solymar,Cogswell (child)  
Results at BTPS

**Technician:**\_\_\_\_\_ **Physician:**\_\_\_\_\_

The first line gives the MicroDL version number and type of report (a summary report of all data stored is also available – see Printing Stored Results).

Subsequent lines give the test number (the tests are numbered sequentially as they are stored), date and time of test, and patient details.

This is followed by the test results presented with the percentage of the predicted value, the predicted value, the predicted minimum and maximum, and the unit of measurement.


The expiratory Flow/Volume and Volume/Time curves are also shown overlaid with the predicted curves and FEV1 marker on the Flow/Volume curve.

To obtain subsequent copies of the report press 'enter'. When all the required printouts have been obtained press 'off' and the unit will turn off.

If the unit is turned off after performing a manoeuvre, and then used again within five minutes, then only the best result within the period will be saved.

## Printing Stored Records

There are five printing options for the stored results. To assist the selection of the required option a system of icons will appear on the display when the unit is set to printing mode:

If one page appears, then a single record has been selected.

If two pages appear, then only records since the last printout will be printed.

If three pages appear, then all stored records will be printed.

If the front page shows only lines of text then a summary report has been selected.

If the Flow/Volume symbol appears, then a spirometry report has been selected.

Depress the up key, turn the unit on and release the up key after a few seconds. The display will indicate the type of report to be printed and the number of pages that will be used.


Use the up key to select the type of printout required.

The following examples assume that the unit is configured in single patient mode, and that symptom scores are recorded.


### **Option 1 – summary report of all records**

After turning the MicroDL on with the up key depressed the following will be displayed:

This display is indicating that a summary report of all stored results has been selected and that the printout will use two sheets of paper.

If this option is required then press 'enter' and the display will change to:

**The following printout of spirometry results will be obtained:**

MicroDL      SUMMARY REPORT

Test	ID	Date	Time	FEV1	FVC	PEF
1	6	10/09/97	10:56	2.97(86%)	3.25(78%)	7.49(87%)
2	6	11/09/97	11:32	2.74(79%)	3.12(75%)	6.45(75%)
3	6	12/09/97	09:58	3.14(91%)	3.53(85%)	8.37(97%)
∨	∨	∨	∨	∨	∨	∨
∨	∨	∨	∨	∨	∨	∨
∨	∨	∨	∨	∨	∨	∨

Values in brackets are percentage of best value.

**Technician:** \_\_\_\_\_ **Physician:** \_\_\_\_\_

**This is followed by the symptom score entries, where applicable:**

MicroDL V1.15      EVENTS REPORT


#	Date	Time	Event	Value
1	10/09/97	10:58	COUGH	8
2	11/09/97	11:34	CAT	
∨	∨	∨	∨	∨
∨	∨	∨	∨	∨
∨	∨	∨	∨	∨

The summary report gives the date, time, FEV1, FVC, and PEF together with the percentage of best value in brackets. The tests are sequentially numbered and this can be used to select a particular test for more detailed investigation by selecting a full report using printing option 2.

Please note that symptom score entry is not possible in multiple patient modes. At any time during printing pressing 'off' will terminate the printing. Should the printer fail to operate correctly then turn off both the printer and the MicroDL unit, check all connections, and try again.


## Option 2 – full spirometry report of one record

After turning the MicroDL on with the up key depressed, use the up key again to obtain printing option 2. The display will show:

This option is used to select a full report, including Flow/Volume and Volume/Time curves, on one stored result. The '1' indicates that one sheet of paper will be used.

Press 'enter' if this option is required and the display will change to:

The tests are numbered sequentially as they are stored and the number displayed represents the last test recorded.

Use the down key to select the required test and then press 'enter'. The display will change to:

The following printout of spirometry results will be obtained:

MicroDL    FULL REPORT


Patient Name:\_\_\_\_\_

Test: 1        Date: 12/11/97    Time:    0802

ID: 6        Sex: Male        Age: 19

Height: 190 cm        Ethnic Correction: 0%

Parameter	Result	%Best	Best	%Pred	Pred	Min	Max	Units
FEV1	3.71	78	4.76	74	4.95	4.11	5.79	L
FVC	4.61	82	5.62	77	5.95	4.95	6.95	L
PEF	9.67	92	10.65	90	10.74	8.75	12.73	L/S
FEV1%FVC	80	97	82	97	82	70	94	%
F50	4.15			68	6.07	3.90	8.24	L/S
F25	1.36			45	2.96	1.68	4.24	L/S
MEF	5.45			102	5.31	3.60	7.20	L/S
FET	4.66							S

Normal Values: ECCS (adult); Zapletal,Solymar,Cogswell (child)

Results at BTPS

Technician:\_\_\_\_\_ Physician:\_\_\_\_\_


The first line gives the MicroDL version number and type of report. Subsequent lines give the test number, date and time of test, and patient details. This is followed by the test results presented with the percentage of personal best value, the personal best value, the percentage of predicted value, the predicted value, the predicted minimum and maximum, and the unit of measurement.

The expiratory Flow/Volume and Volume/Time curves are also shown overlaid with the predicted curves and FEV1 marker on the Flow/Volume curve. Please note that personal best values cannot be entered in multiple patient modes.

To obtain subsequent copies of the report press 'enter'. When all the required printouts have been obtained press 'off' and the unit will turn off.


### **Option 3 – full spirometry report of all records**

After turning the MicroDL on with the up key depressed, use the up key twice to obtain printing option 3. The display will show:

This option is used to print full spirometry reports, including Flow/Volume and Volume/Time curves, on all the stored results.

The '65' indicates that sixty-five sheets of paper will be used, as there are currently sixty-five tests stored. Press 'enter' if this option is required. The display will show:


The printouts obtained will be identical to the single report produced with option 2.

At any time during printing pressing 'off' will terminate the printing. Should the printer fail to operate correctly then turn off both the printer and the MicroDL unit, check all connections, and try again.

#### **Option 4 – summary report of all records since last printout.**

After turning the MicroDL on with the up key depressed, use the up key three times to obtain printing option 4.

The display will show the following:

This option is used to print a summary report of all the results that have been recorded since the last summary report was printed.

Press 'enter' if this option is required. The display will show:


The printouts obtained will be identical to the summary report produced with option 1, but with fewer entries if a summary report has been previously obtained.

### **Option 5 – full spirometry report of all records since last printout.**


After turning the MicroDL on with the up key depressed, use the up key four times to obtain printing option 4.

The display will show:

This option is used to print full spirometry reports of all the results that have been recorded since the last time that printing option 3 was used.

Press 'enter' if this option is required. The display will show:

The printouts obtained will be identical to the single report produced with option 3.

## **PC Based Spirometry**

The MicroDL can be connected to a PC and interfaced to CareFusion spirometry data management software SPCS. SPCS enables you to record, manage, and display spirometry information. Its main features are:

A fully featured spirometry database.

An interface, which enables the unit to be controlled from the PC and provides live display of Flow/Volume, and Volume/Time graphs together with a child incentive display.


Reporting facility including a full set of parameters, predicted values, graphs and interpretations.

Trend analysis of spirometry data.

## **Configuring the MicroDL with SPCS**

The MicroDL can be set up using SPCS. First, ensure that the MicroDL is connected to the PC using the cable supplied. Connect the nine way 'D' connector to the serial port of the PC and the 3 pole jack plug to the socket in the right hand side of the MicroDL.

To communicate with the PC the MicroDL must be turned on in PC mode. Depress the 'enter' key, turn the unit on, and release the 'enter' key when the display shows:

Start SPCS on the PC. Ensure that SPCS is set up to use the COM port on the computer that the MicroDL is connected to (refer to online SPCS help). A patient record should then be set up containing the details to be downloaded to the MicroDL (or an existing patient selected). These details can then be used to initialise the unit. Please refer to the online help included with SPCS for the exact details for performing these operations.

NOTE: Keep the PC out of reach of the patient when connected to the MicroDL. The MicroDL should only be connected to a computer that is manufactured in accordance with EN60950 1992/1993 'Safety of Information Technology Equipment including Electrical Business Equipment'.

## Uploading Results

To upload the results stored in the MicroDL, connect the unit to the PC as described above and start SPCS. When the initial screen is displayed, click on the 'upload' button. The data contained in the unit will then be uploaded into the SPCS database (refer to online help for further information).

Please note that the patient details for any data uploaded into SPCS from a unit that has been set-up by hand, should be entered into the SPCS database. This will allow the uploaded spirometry records to be displayed and trended.

The patient IDs used in the unit should match those in the SPCS database.

NOTE: Keep the PC out of reach of the patient when connected to the MicroDL.

To connect the unit to the PC; connect the cable to one of the PC's serial communications ports and turn the unit on with the 'enter' key held down. The unit display should show 'PC' flashing. Now run SPCS on your PC (for more information on SPCS, refer to online help).

## Battery Management

There are two sources of internal battery power in the MicroDL.

The recorded data is stored in a battery-backed memory, supplied by a Lithium power source with an expected life of at least 10 years. This battery also supplies the timekeeping clock.

The microcomputer circuitry is supplied by a 9 volt, alkaline PP3 battery. The battery life is approximately 24 hours of continuous use.

The battery voltage is monitored continuously and the following warning message is displayed when the battery starts to become exhausted:

The warning will flash for a few seconds and the unit will turn itself off. Replace the battery by sliding off the battery cover at the back of the unit, unclip the exhausted battery and replace with a new battery of the correct type.

Note: if the battery is not to be used for long periods, the battery should be removed to prevent damage to the instrument by possible leakage.

## Calibration Check

The MicroDL is calibrated to read in litres at body temperature, barometric pressure, and saturated with water vapour (BTPS) for FEV1 and FVC. The units of PEF are litres per second BTPS or litres

per minute BTPS depending upon the configuration in the PC software.

The calibration should remain stable indefinitely, unless the transducer is physically damaged, and the unit should not require re-calibration.

However, as a check on the correct functioning of the unit we do recommend that the calibration is checked periodically.

To check the calibration of the unit a large syringe should be used.

This should have a 30mm outlet and be of at least three litres in capacity and preferably six litres.

A 3 Litre Syringe is available from CareFusion. The MicroDL should be connected to this syringe with the minimum of adapters and the syringe volume should be injected into the transducer evenly, without pausing.

If, after checking, the unit appears to require re-calibration then please refer to your supplier for service.

## **Looking after your MicroDL Spirometer**

Please observe the following precautions:

- \* Avoid exposing the MicroDL to direct sunlight.
- \* Avoid operating the spirometer in dusty conditions or near to heating appliances or radiators.
- \* Do not keep the spirometer in a damp place or expose it to extremes of temperature.
- \* Do not direct the transducer holder towards a strong light source whilst operating the spirometer.

## Cleaning the Transducer

The transducer requires no routine maintenance or servicing. However, if you wish to sterilise or clean the transducer it may be removed by means of the following procedure:

1. Remove the transducer by gently pulling from the main body with a twisting action.
2. The transducer may now be immersed in warm soapy water for routine cleaning or immersed in cold sterilising solutions e.g. Perasafe, for a period not exceeding 10 minutes (Alcohol and chloride solutions should be avoided). After cleaning/sterilising, the transducer should be rinsed in distilled water and dried.  
Perasafe is available from CareFusion in convenient 81g plastic bottles.
3. Re-assemble the transducer into the MicroDL.

## Caution

- Pulmonary filters are single patient use. If used on more than one patient there is a risk of cross-infection. Repeat use may increase air resistance and lead to an incorrect measurement.”
- Mouthpieces are single patient use. If used on more than one patient there is a risk of cross-infection. Repeat use may degrade materials and lead to an incorrect measurement.”

## Symbols

Type B device

## Environment

This instrument complies with directive EN60601-1-2 electromagnetic compatibility but can be affected by cellular phones and by electromagnetic interference exceeding levels specified in EN 50082-1:1992

## Servicing

A full service manual including circuit diagrams and parts list is available upon request.

For US Customer Care: Toll Free 1-800-231-2466, phone 1-714-283-2228

## Specification of the MicroDL

### Displayed Measurements:

Forced Vital Capacity (FVC)  
Forced Expired Volume in 1 second (FEV1)  
Peak Expiratory Flow Rate (PEF),  
Forced Expiratory Ratio (FEV1 FVC)

### Recorded Measurements:

Forced Vital Capacity (FVC)  
Forced Expired Volume in 1 second (FEV1)  
Peak Expiratory Flow Rate (PEF),  
Forced Expiratory Ratio (FEV1%FVC)  
Flow at 50% of volume remaining (F50)  
Flow at 25% of volume remaining (F25)

Mid Expiratory Flow (MEF)  
Forced Expiratory Time (FET)

**Storage Capacity:** Typically > 1 year of data recording

**Display:** Custom LCD

**Transducer Type:** Uni-Directional Digital Volume

**Accuracy:** To ATS recommendations  
Standardisation of spirometry

**Data Entry:** 4 event markers or symptom scores on  
1 to 10 scale

**Power supply:** 9 volt PP3 Alkaline cell

**Backup battery life:** >10 years

**Main Battery Life:** Approximately 24 hours of continuous  
use

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## Customer contact information

### **UK Customers only**

For all Sales Order processing for products, training and spare parts, Service and Technical Support enquiries, please contact the following:

CareFusion UK 232 Ltd  
UK Customer Service & Support  
The Crescent  
Jays Close  
Basingstoke  
RG22 4BS

Customer Service Sales Enquiries:

Telephone: 01256 388550

Email: [micro.uksales@carefusion.com](mailto:micro.uksales@carefusion.com)

Factory Repair and Administration Enquiries:

Telephone: 01256 388552

Email: [micro.ukservice@carefusion.com](mailto:micro.ukservice@carefusion.com)

Technical Support Enquiries:

Telephone: 01256 388551

Email: [support.rt.eu@carefusion.com](mailto:support.rt.eu@carefusion.com)

### **International customers only**

For all Sales Order processing for products and Spare parts, Service and Technical Support enquiries, please contact the following:

Carefusion Germany 234 GmbH  
Customer Service & Support International  
Leibnizstrasse 7  
97204 Hoechberg  
Germany

Customer Service Sales Enquiries:

Telephone: 0049 931 4972 670

Email: [micro.internationalsales@carefusion.com](mailto:micro.internationalsales@carefusion.com)

Factory Repair and Administration Enquiries:

Telephone: 0049 931 4972 867

Email: [support.admin.eu@carefusion.com](mailto:support.admin.eu@carefusion.com)

Technical Support Enquiries:

Telephone: 0049 931 4972 127

Email: [support.rt.eu@carefusion.com](mailto:support.rt.eu@carefusion.com)

CareFusion UK 232 Ltd.,  
Quayside  
Chatham Maritime  
Kent ME4 4QY  
U.K.

**CE**  
0086

JDE no. 36-MAN1283  
Drg no. 076-06WW  
Issue 1.0  
December 2009  
© CareFusion 2009